

Introduction

- Food samples like sunflower oil need to be analyzed for unwanted compounds
- Polycyclic aromatic hydrocarbons are cancerogenic compounds from incomplete combustion of organic matter
- Sample needs to be cleaned-up in a way that enables the injection to a suiting instrument like a gas chromatography mass spectrometer (GC-MS).
- "Clean-up" especially from edible oil and fat samples is complex and labor-intensive
- For a sample like sunflower oil, all lipids need to be removed before GC-MS injection.
- However, clean-up with established methods is labor intensive

PAH in oil samples?



The greater concepts



The Solution: Automation with μ -SPE

- Automation is one of the strategies to increase efficiency of the SPE clean-up
- Automated SPE methods using mini-tubes, so called μ -SPE have already been successfully implemented for clean-up of QuEChERS extracts used for analysis of PAH and other analyte classes such as pesticides and PCB's.
- The miniaturization of the clean-up step to a microliter scale solid phase extraction (μ SPE).The PAH fraction is eluted only in a small volume of few 100 μ L for direct injection into GC-MS
- Extraction and evaporation with a potential loss of compounds is avoided

Objective

The objective was to establish a μ -SPE clean-up procedure for PAH analysis in sunflower oil. We were using two μ -SPE cartridges with Florisil and C18/Z-Sep. We were using the greater concepts of Green Chemistry and Automation.

μ -SPE: The general concept




Figure 1: PAL RTC μ SPE Clean-up Configuration

μ -SPE: The developed concept




Figure 2: Concept of the Automated μ SPE Clean-up Steps

Chromatogram and detailed procedure




Figure 3 : GC-MS total ion current (TIC) chromatogram of a sunflower oil eluent after using μ SPE clean-up.

Sequence	Description	Tool
1	Change Tool to 1ml syringe	1ml syringe
2	Wash syringe two times with 1000 μ l acetonitrile	1ml syringe
3	Load 240 μ l of acetone to cartridge 1 (florisil) for conditioning	1ml syringe
4	Drying with air, 1000 μ l, cartridge 1	1ml syringe
5	Change Tool to 10 μ l syringe (Injection Tool)	1ml syringe /Injection Tool
6	Wash syringe two times with 10 μ l acetonitrile	Injection Tool
7	Load 4 μ l of oil sample to the syringe (2 fill strokes prior to avoid bubbles)	Injection Tool
8	Load extracted fat sample to cartridge 1 (florisil)	Injection Tool
9	Wash injection tool two times with 10 μ l acetonitrile	Injection Tool
10	Change Tool to 1ml syringe	1ml syringe
11	Wash syringe two times with 1000 μ l acetonitrile	1ml syringe
12	Elute cartridge to vial 1 with 360 μ l acetonitrile	1ml syringe
13	Wash syringe two times with 1000 μ l acetonitrile	1ml syringe
14	Load 240 μ l of Acetone Condition of cartridge 2	1ml syringe
15	Drying cartridge 2 with 1000 μ l air	1ml syringe
17	Load/Eluate 260 μ l of eluate 1 to cartridge 2 above vial 2	1ml syringe
18	Wash syringe two times with 500 μ l acetonitrile	1ml syringe

Figure 4: The Automated μ SPE Clean-up Steps

Validation with 24 PAH + 22 IS

- Analytes: 24 PAH analytes
- Internal Standards: 22 PAH
- 5 replicates were used
- Sample: Sunflower oil
- Our LOQ: 1090 μ g/kg
- Maximum level (sum of 4 PAH): 10 μ g/kg




Figure 5: Direct PAL clean-up and analysis

<https://onlinelibrary.wiley.com/doi/10.1002/jssc.202000720>

Total ion chromatogram after clean-up with μ SPE

Table 1: Results from recovery analysis of sunflower oil spiked at two levels (326 and 3260 μ g/kg) in five replicates.

Analyte	Recovery (RSD)
5-Methylchrysene	98(3)
Acenaphthylene	101(4)
Acenaphthene	107(8)
Anthracene + Phenanthrene	99(5)
Benz[a]anthracene + Chrysene	90(6)
Benzo[a]pyrene	104(13)
Benzo[c]fluorene	93(5)
Benzo[g,h,i]perylene	73(22)
Benzo[b+j+k]fluoranthene	115(2)
Cyclopenta[cd]pyrene	105(3)
Dibenz[a,h]anthracene	96(8)
Dibenzo[a,e]pyrene	72(19)
Dibenzo[a,L]pyrene	69(17)
Fluoranthene	107(5)
Fluorene	102(5)
Indeno[1,2,3-cd]pyrene	53(14)
Naphthalene	97(8)
Pyrene	97(5)

Conclusion

- New micro-SPE method developed
- Vegetable oil sample could be directly cleaned up
- PAH could be analyzed with good recoveries
- Improvements on sensitivity can be used for decreasing LOQ
- Larger volume injections can be applied for increased sensitivity
- Fully green analytical method, evaporation steps avoided

Presented at

LAPRW Conference 2021, May 18-20, Online-event.

References

- Valentina Bertoz, Giorgia Purcaro, et al., A Review on the Occurrence and Analytical Determination of PAHs in Olive Oils., Foods 2021, 10(2), 324; <https://doi.org/10.3390/foods10020324>.
- Bruce D. Morris and Richard B. Schriener, Development of an Automated Column Solid-Phase Extraction Cleanup of QuEChERS Extracts, Using a Zirconia-Based Sorbent, for Pesticide Residue Analyses by LC-MS/MS, J. Agric. Food Chem. 2015, 63, 5107-5119, DOI: 10.1021/ef505539e.
- Steven J. Lehotay, Lijun Han, Yelena Sapozhnikova, Automated Mini-Column Solid-Phase Extraction Cleanup for High-Throughput Analysis of Chemical Contaminants in Foods by Low-Pressure Gas Chromatography-Tandem Mass Spectrometry, published with open access at Springerlink.com, Chromatographia (2016) DOI 10.1007/s10337-016-3116-y.
- Eyring P, Preiswerk T, Frandsen HL, Duedahl-Olesen L (2020) Automated micro-solid-phase extraction clean-up of polycyclic aromatic hydrocarbons in food oils for analysis by gas chromatography – orbital ion trap mass spectrometry. J Sep Sci 1-9. <https://doi.org/10.1002/jssc.202000720>