



The New PAL Micro-SPE Cartridge for Automated Pesticides Extraction and Clean-up

Micro-SPE for Green Analytical Chemistry


PAL SYSTEM
Ingenious sample handling

Solid Phase Extraction

A quick look back and where do we go today

- ‘The separating funnel is a museum piece’
- Modern SPE originated in 1974
by Reginald Adams, Thomas Good, and Michael Telepchak
First dispersive (dSPE)
Later cartridge formats
- Much simpler for the lab
 - Less sorbent material
 - Less solvent
 - Faster
 - More concentrated analytes
 - More selective
 - Compatible with GC-MS and LC-MS

Exactly the same is true again for μ SPE




What is Micro-SPE (μ SPE)

Compare to the classical cartridge SPE

Classical SPE


- Limited selectivity
 - High sample and solvent volumes
 - Requires evaporation with N_2
 - End volume $\gg 100 \mu\text{L}$ in vial
- Vacuum operated
- Drying before elution
- Manual operation
 - Time consuming
 - Low sample throughput
 - Batch processing
- No QA/QC
 - As of manual operation



μ SPE


- High selectivity
 - Sharp elution peak profile,
 - Compares to LC separation
 - No dilution, no concentration
 - Final eluates $< 100 \mu\text{L}$ (or online)
- Positive pressure w liquid syringe
 - Very low solvent use
- No drying step
- Walk away automation
 - Fast
 - Works on chromatographic timescale
 - High productivity
- Traceable
 - Processing well documented

*Works like a
mini LC column !*




Both available for μ SPE

- “Enrichment” mode (the “classical” procedure)
aka Load-Wash-Elute mode
Analytes ● retained
Matrix ▲■ washed away
Analytes ● eluted by solvent change
e.g. with C18, or SAX material for Glyphosate, AMPA in EURL Almeria
- “Scavenging” mode
Matrix ▲■ retained
kept on cartridge
Analytes ● elute with extract
e.g. for ‘QuEChERS’ with acetonitrile, or ‘SweEt’ with ethyl acetate
 - Pesticides in Hill Labs, USDA Lehotay, Zurich Kanton Lab, EURL Almeria ...)
 - C18 material for veterinary drugs analysis



Inside the μ SPE Cartridge

How does the μ SPE Cartridge work




- The μ SPE cartridges offer combinations of sorbent materials - as used for the QuEChERS clean-up, customized and proprietary sorbents are available, just filter materials of different pore sizes, e.g. for LC and IC applications.

PAL System μ SPE – Automated Cartridge Assembly

Automated cartridge filling and assembly with built in QC



• PAL System μ SPE cartridges components




Automated QC reporting and checks

- Production of sorbent mixture (dry storage)
 - Testing on the dosed amount of powder (weight)
 - Check for correct packing density of the powder (height)
 - Lot number assignment
 - Packing and labelling of cartridges from a produced lot
 - Dry storage of packed cartridges
-
- Lot tracking of production
 - Sorbent amount (weight)
 - Sorbent packing (height)
-
- Certificate of Conformity for each lot
 - Filling Amount: Target weight $\pm 10\%$
 - Overall cartridge height: Target height $\pm 1.0\text{mm}$



Production Control



		Statistics	
S	Mean Value and Deviation	Average CF	15.13 mg
		Std. Dev.	1.00 mg
G1	Percentile Quality level limits	Perc. 0.05	14.70 mg
		Perc. 0.95	15.47 mg
E	Limits- / Extreme values	Min. Value	-0.10 mg
		Max. Value	25.00 mg
F	Out of Quality Reasons	Out of Tolerance	58 pcs
		Missing Filter	22 pcs
		IBS Height Out of Tol.	1 pcs

PAL System μ SPE - QC and Reporting

Automated QC reporting and checks

- Production of sorbent mixture (dry storage)
 - Testing on the dosed amount of powder (weight)
 - Check for correct packing density of the powder (height)
 - Lot number assignment
 - Packing and labelling of cartridges from a produced lot
 - Dry storage of packed cartridges
-
- Lot tracking of production
 - Sorbent amount (weight)
 - Sorbent packing (height)



Product
PAL SYSTEM Certificate Of Conformity

Part number	uSPE-GCQuE1-45-V	
Description	uSPE Kartusche GC QuEChERS sorbent mix 1.45mg	
Cartridge type	01-05B	
Content	MgSO4	20 mg
	PSA	12 mg
	C18-EC	12 mg
	GCB	1 mg

The μ SPE cartridges are manufactured and tested to comply with the following specifications:

	Specification	QC
Target filling amount	45 mg	Pass
Filling tolerance	$\pm 10\%$	Pass
Target overall height	36.4 mm	Pass
Overall height tolerance	± 0.5 mm	Pass


Christian Franz
Manager QA
This document has been generated electronically and is valid without signature.

Revision 1
TPA-000074 COC_USPE


15.13 mg
1.00 mg
14.70 mg
15.47 mg
-0.10 mg

How Does μ SPE Work on a PAL System

One Trayholder is installed on a PAL System (upgradeable)



- ① Raw Extracts
e.g. QuEChERS, SweEt, or other ...
- ② Elution - Cleaned Extracts
for injection GC-MS, LC-MS
- ③ μ SPE cartridges
for opt. conditioning



Cartridge transport on syringe needle
for elution and return/waste bin

Standard QuEChERS Protocol*

10 g of sample

Shake 1 min

Salt out

Shake 1 min

Centrifuge

Dispersive SPE

Shake 1 min

Centrifuge

Analysis

Use 50 mL tubes
Add 10 mL acidified acetonitrile
Add ISTD (Triphenylphosphate)

4 g MgSO₄ anh.,
1 g NaCl,
1 g Citrate buffer (CEN 15662)
1 g Acetate buffer (AOAC 2007.01)

- ➔ Freeze step for fatty samples, or GLP clean-up
- ➔ Direct LCMS analysis of polar pesticides

1 mL aliquot transferred into 10 mL tubes,
add PSA, MgSO₄, C18, GCB, Chlorofiltr, ...
as required by the food commodity

LLE Extraction



dSPE Clean-up



* QuEChERS - Mini-Multiresidue Method for the Analysis of Pesticides, M. Anastassiades, 2003

QuEChERS Protocol* with μ SPE Clean-up

10 g of sample

Shake 1 min

Salt out

Shake 1 min

Centrifuge

~~Dispersive SPE~~

~~Shake 1 min~~

~~Centrifuge~~



Use 50 mL tubes
Add 10 mL acidified acetonitrile
Add ISTD (Triphenylphosphate)

➔ Add buffer kit, AOAC 2007.1/EN15662

➔ No freeze - No GPC

➔ Take 1 mL for μ SPE clean-up

Black and green tea samples



Chlorpyrifos 1 ppb
RT 13.45 min
m/z 313.90 > 285.90
m/z 313.90 > 193.90

Parathion 1 ppb
RT 13.59 min
m/z 291.10 > 137.00
m/z 291.10 > 81.00

GC-MS and LC-MS Injection

Pesticides Clean-up by Hill Laboratories, Hamilton, NZ

More than 10 Years Experience for > 900 Matrices, >1000 Samples/wk

- Automated, compared to original QuEChERS
- Reduced manual labor in sample prep
- Wider range of samples
high lipid content, incl. Avocados, with difficult matrices like dried herbs and spices
- Off-line clean-up on PAL Systems
only 1/5 of analysis times, serves several GC-MS and LC-MS systems.



μSPE GC-MS clean-up		
Sorbent	Amount	Percentage
PSA	12 mg	27%
C18	12 mg	27%
GCB	1 mg	2%
MgSO ₄	20 mg	44%
Total	45 mg	100%

μSPE LC-MS clean-up		
Sorbent	Amount	Percentage
Z-Sep	8 mg	27%
C18	21 mg	70%
GCB	1 mg	3%
Total	30 mg	100%

JOURNAL OF AGRICULTURAL AND FOOD CHEMISTRY Article
pubs.acs.org/JAFAC

Development of an Automated Column Solid-Phase Extraction Cleanup of QuEChERS Extracts, Using a Zirconia-Based Sorbent, for Pesticide Residue Analyses by LC-MS/MS

Bruce D. Morris* and Richard B. Schriener

Food and Bioanalytical Division, R. J. Hill Laboratories, Private Bag 3205, Hamilton East, New Zealand

Supporting Information

ABSTRACT: A new, automated, high-throughput, mini-column solid-phase extraction (c-SPE) cleanup method for QuEChERS extracts was developed, using a robotic X-Y-Z instrument autosampler, for analysis of pesticide residues in fruits and vegetables by LC-MS/MS. Removal of avocado matrix and recoveries of 263 pesticides and metabolites were studied, using various stationary phase mixtures, including zirconia-based sorbents, and elution with acetonitrile. These experiments allowed selection of a sorbent mixture consisting of zirconia, C₁₈, and carbon-coated silica, that effectively retained avocado matrix but also retained 53 pesticides with <70% recoveries. Addition of MeOH to the elution solvent improved pesticide recoveries from zirconia, as did citrate ions in CEN QuEChERS extracts. Finally, formate buffer in acetonitrile/MeOH (1:1) was required to give >70% recoveries of all 263 pesticides. Analysis of avocado extracts by LC-Q-Orbitrap-MS showed that the method developed was removing >90% of di- and triacylglycerols. The method was validated for 269 pesticides (including homologues and metabolites) in avocado and citrus. Spike recoveries were within 70–120% and 20% RSD for 243 of these analytes in avocado and 254 in citrus, when calibrated against solvent-only standards, indicating effective matrix removal and minimal electrospray ionization suppression.

KEYWORDS: QuEChERS, SPE, LC-MS/MS, ITSP, Z-Sep, zirconia, pesticide, multiresidue, avocado, citrus

INTRODUCTION

The "quick, easy, cheap, effective, rugged, and safe" (QuEChERS) method for the analysis of multiclass pesticide residues in fruits and vegetables introduced the use of dispersive solid-phase extraction (d-SPE) cleanup, to reduce the amounts of coextracted matrix in extracts, before instrumental analysis, using a mixture of MgSO₄ and PSA sorbents, with C₁₈ and graphitized carbon black (GCB) added if required to improve removal of nonpolar matrix and diolophyll.^{1–4} In the original QuEChERS method, d-SPE was used instead of column SPE (c-SPE) to provide a quicker and cheaper cleanup.¹ Recently the zirconia-based sorbent HybridSPE, in well-plates or columns, has been utilized for the removal of phospholipids from plasma^{5–7} and eggs.⁸ The zirconia materials Z-Sep and Z-Sep+ have been evaluated for d-SPE cleanup of QuEChERS extracts for analysis of environmental pollutants and pesticides in fish and shrimp^{9–11} and pesticides from oily fruits or vegetable oils,^{12–16} due to their abilities to remove the lipophilic matrix. However, in our experience, used routinely, Z-Sep d-SPE can result in the transfer of solid phase into analysis vials and subsequently into the HPLC, building up over time to cause retention of some analytes and poor peak shapes or carry-over. Consequently, we investigated the development of an automated c-SPE cleanup, based on zirconia-coated silica, using Instrument Top Sample Preparation (ITSP) mini-cartridges, on a robotic X-Y-Z instrument autosampler. This could be as quick and cheap as d-SPE, as many instruments are already equipped with robotic autosamplers; however, it could also give the improved matrix removal that is possible with c-SPE¹⁷ and avoid zirconia transfer to the LC-MS/MS.

Avocado extracts were selected as a matrix with high oil content,^{18,19} and experiments were carried out to evaluate the weight of matrix removed after acetonitrile (MeCN) elution through ITSP c-SPE cartridges with six different stationary phases. Recoveries of 263 pesticides and metabolites spiked on avocado were determined through five of these sorbents and, along with matrix weight-removal results, allowed selection of a Z-Sep/C₁₈/CarbonX mixture for further method development. Investigation of the effect of different elution solvents (MeCN, MeCN/MeOH (1:1), MeOH, coextracted citrate in a CEN (European Committee for Standardization method),²⁰) QuEChERS extract, and formate buffer at three concentrations in MeCN/MeOH (1:1), on pesticide recoveries through Z-Sep/C₁₈/CarbonX, resulted in a method using elution of CEN QuEChERS extracts with 100 mM formate buffer in MeCN/MeOH (1:1). To the best of our knowledge, this study is the first to use ITSP mini-cartridges for cleanup of QuEChERS extracts and zirconia solid phase in an SPE column, rather than used despectively, for pesticide residue analysis. Removal of avocado di- and triacylglycerols by Z-Sep, monitored by LC-Q-Orbitrap-MS, is also presented. The method was validated for the analysis of 269 pesticides, including homologues and metabolites, in avocado and citrus, to give spike recovery and reproducibility data.


Special Issue: 51st North American Chemical Residue Workshop
Received: November 17, 2014
Revised: February 7, 2015
Accepted: February 9, 2015
Published: February 23, 2015

ACS Publications © 2015 American Chemical Society 5107
DOI: 10.1021/acs.jafc.5b0119
J. Agric. Food Chem. 2015, 63, 5107–5119

B.Morris, R.Schriener (2015) J. Agri. Food. Chem. 63, 5107-19.

Automated μ SPE Clean-up for GC-MS - Only 8 min *

1. Take a 1 mL syringe
2. Wash the 1 mL syringe with MeCN
3. Take 500 μ L raw extract from rack 1 into 1 mL syringe
4. Get the μ SPE cartridge from rack 3 with the syringe needle
5. Move the cartridge to the elution rack 2
6. Push the raw extract through the μ SPE cartridge at 5 μ L/s
7. Discard μ SPE cartridge into waste beaker
8. Wash the 1 mL syringe with MeCN/MeOH/water (vol 1/1/1)
9. Wash the 1 mL syringe with MeCN
10. Switch to 100 μ L syringe and wash with MeCN
11. Add 25 μ L AP + QC solutions to the collection vial in rack 2
12. Wash the 100 μ L syringe with MeCN/MeOH/water (vol 1/1/1)
13. Wash the 100 μ L syringe with MeCN
14. Switch to 10 μ L GC injection syringe
15. Wash the 10 μ L syringe with MeCN
16. Aspirate the cleaned extract from the elution vial in rack 2
17. Inject 1 μ L of extract to GC-MS/MS
18. Wash the 10 μ L syringe with MeCN




* as of Steven J. Lehotay, Lijun Han, Yelena Sapozhnikova (2016) and Nicolas Michlig, Steven J. Lehotay (2022)

Highly Polar Pesticides in Complex Matrices (QuPPE)

Glyphosate, AMPA, Glufosinate, ... EURL for Pesticides in Fruit and Vegetables, Almeria, Spain

- Matrix honey, pollen, coffee beans
- Acidified methanol extraction
- Automated μ SPE clean-up
50 mg SAX (strong ion exchange)
- Clean-up procedure
1000 μ L methanolic raw extract
Load at 5 μ L/s
Matrix washed with 600 μ L methanol
Analytes elute with 400 μ L methanol/HCl (9:1)
- Inject 10 μ L to LC-MS/MS
- Cost saving: 500 mg \rightarrow 50 mg SAX material
- Time saving: manual \rightarrow automated 10:1
- Analytical: Improved recoveries up from avg. 70 \rightarrow 86%




Extracted LC-MS/MS ion chromatograms of standards of AMPA, glufosinate, and glyphosate at 0.010 mg/kg spiked to honey and pollen matrix without (left) and after SAX clean-up (right)

Jesus, F., A. R. Garcia, et al. 2023. "Determination of Highly Polar Anionic Pesticides in Beehive Products by Hydrophilic Interaction Liquid Chromatography Coupled to Mass Spectrometry." *Anal. Bioanal. Chem.* <https://doi.org/10.1007/s00216-023-04946-7>.

Ethylacetate Extraction of Pesticides from Foods (aka SweEt)

Big time savings and reduced manual effort for high fat matrices – Cantonal Lab. Zurich, CH

- EtOAc extracts wider range of polar pesticides
- But, also extract high amounts of matrix
- GPC or extract freezing was used as clean-up
- Clean-up using μ SPE
 - 45 mg of PSA, C18, GCB, $MgSO_4$
 - Load 200 μ L raw EtOAc extract, 2 μ L/s
 - Blow-out 1 mL air
- Injection 3 μ L cleaned extract to GC-MS/MS
- Significant improvements:
 - Difficult and fatty samples e.g. dried spices, egg, avocado, or liver are successfully processed,
 - Captan and Folpet with good recoveries.
 - One μ SPE cartridge for all food matrices.
 - No time-consuming freeze-out or GPC required anymore.



Captan and folpet in raspberry samples after μ SPE clean-up.

Schürmann, A., C. Crüzer, et al. 2023 "Automated Micro-Solid-Phase Extraction Clean-up and Gas Chromatography-Tandem Mass Spectrometry Analysis of Pesticides in Foods Extracted with Ethyl Acetate." Anal. Bioanal. Chem. 416 (3): 689–700. <https://doi.org/10.1007/s00216-023-05027-5>.



Veterinary Drug Screening by online μ SPE – LC-MS/MS

From the Veterinary Diagnostic Laboratory, Iowa State University, USA.

- Veterinary drugs are legally controlled in large number of samples in difficult matrix.
- QuEChERS extraction (LLE with MeCN) from 5 g muscle, or kidney
- Automated μ SPE clean-up on Thermo Scientific™ TriPlus™ RSH sampler, μ SPE cartridge with 15 mg of endcapped C18. 300 μ L of supernatant at 2 μ L/s, the eluate diluted 3+1 with mobile phase, injected into a 2 μ L loop on the injection valve. Clean-up takes only 8.5 min to complete
- Cost saving on C18: 500 mg \rightarrow 15 mg (30x less)
- Time saving: 80 min/15 samples \rightarrow Zero, online prep
- No additional consumables




103 veterinary drugs at 50 ng/g in bovine kidney extract. Total cycle time 15 min.



Automated QuEChERS Extraction and μ SPE Clean-up

Application for homogeneous samples like beverages

- Why not automate all steps? The QuEChERSs extraction with SPE clean-up?
- Homogenous samples
 - Do not need much manual treatment
 - Can be pipetted into 2 mL vials (also automatically!)
- QuEChERS extraction is automated in 2 mL vials with buffer salts previously added. Then added acidified MeCN, sat. NaCl solution.
- μ SPE clean-up
 - 45 mg of PSA, C18, GCB, MgSO₄
 - Load 250 μ L raw extract, 2 μ L/s
- Injection 3 μ L cleaned extract to GC-MS/MS
- Combines all benefits from QuEChERS, μ SPE clean-up, and prep-ahead automation



PAL RTC System for QuEChERS, μ SPE, and online GC-MS



Tool Park Station

- Pos.1 GC Injection syringe (10 μ L)
- Pos.2 μ SPE Tool (1000 μ L)
- Pos.3 APs/ISTD syringe (25 μ L)





Solvent Module

- Pos.1 Acetonitrile
- Pos.2 NaCl, sat.
- Pos.3 not used

Fast Wash Module

- Pos.1 Water
- Pos.2 Acetonitrile


QuEChERS Extraction and μ SPE Clean-up Workflow

Liquid-Liquid Extraction	> Use Tool	1 Prepare 500 μ L of Orange Juice into 2mL vial	 Orange juice  Orange juice + acetonitrile  + NaCl for phase separation  After μ SPE clean up, for injection
	> Repeat		
	> Get Liquid From Vial	2 Add 600 μ L MeCN into Orange Juice sample	
	> Dispense Liq. Into Vial		
	> Repeat End	Shake	
	> Vortex Vial		
	> Clean Syringe		
	> Get Liquid From Vial	Add 250 μ L sat. NaCl into sample vial, shake	
	> Dispense Liq. Into Vial		
	> Clean Syringe	3 Wait for phase separation from the liquid-liquid extraction	
> Vortex Vial			
> Wait			
μ SPE	> Use Tool	4 Load 250 μ L of the orange extract to the μ SPE clean-up	
	> μ SPE Elute		
	> Clean Syringe		
Inject	> Liquid Injection	Inject 2 μ L of the cleaned extract to GC-MS/MS	
	> Wait		

OCP Pesticides PAL Extracted from Orange Juice


Total ion chromatogram from GC-MS/MS in MRM mode

- At MRL level the pre-spike RSDs were mostly below 10%. Recoveries achieved 70% to 115%, calibration linearity > 0.995.
- LODs range from 1.8 ng/mL to 4.1 ng/mL (n = 8) well below the general MRL at 10 ng/g level.
In the original orange juice from a local supermarket about 1.8 ng/mL Malathion was detected.



Parallel Extraction and μ SPE Clean-up by PAL Prep-ahead

PAL serves sample prep and GC-MS/LC-MS analysis in parallel



- Strong benefits:
- High reproducibility: All samples are treated on the same timeline.
- Saves time: Continuous analysis of samples, improved sample throughput, overnight processing.
- Highly efficient: Increased use of the GC-MS and LC-MS unit, reduces cost/sample.

μSPE replaces the traditional SPE concentration and clean-up procedures

- μSPE is the next step available towards a Greener Analytical Chemistry

- Less solvents
- Less consumables
- Less waste
- Less energy consumption

- μSPE delivers strong analytical advantages

- One clean-up cartridge for all type of samples
- Improved recoveries
- Improved clean-up
- Improved precision

- μSPE reduces cost/sample

- Efficient use of GC-MS and LC-MS by online prep-ahead
- Increased sample throughput
- Walk-away automation
- Less manual workload
- Less repeat measurements
- Faster report out



In Vendor Seminars, Presentations and Posters

Today is Theme Day: Green Analytical Chemistry, Miniaturization and Automation

- **V-5 Gerstel Vendor Seminar**, Wednesday, 18 September 2024 | 10:50 – 11:15
- Steven J. Lehotay, Nicolas Michlig: Comparison of Automated vs. Centrifugal μ-SPE for the Analysis of Pesticides and Other Contaminants in a Variety of Foods.
- **O-16 Steven J. Lehotay**: High Throughput, Miniaturized, Automated, and “Green” Analysis of a Wide Scope of Pesticides and Other Residues in Fatty and Non-fatty Foods.
- **O-17 Amadeo Fernandez-Alba**: Automation and Miniaturization for Enhancing Analytical Methods in Pesticide Residue Evaluation in Food.
- **O-27 Mette Erecius Poulsen**: Evaluation, Optimization, and Validation of Different Micro-SPE Clean-up Cartridges for Pesticide Residuals in Cereals.

And a number of interesting posters:

- **PD-01 Jose Antonio Martinez**, et al. Pesticide Residue Analysis in Spice Samples by Automated Micro-SPE Clean-up
- **PD-18 Ivan Aloisi** et al. High Throughput Automated Off-line Micro-SPE for Clean-up of Complex QuEChERS Extracts
- **PD-49 Ederina Ninga**, et al. Optimization of Customized Micro-SPE Cartridges for Pesticides in Cereals
- **PD-50 Ederina Ninga**, et al. Validation of Pesticides in Insects Using “QuEChERSER”-Micro-SPE Clean-up
- **PD-51 Ederina Ninga**, et al. QuEChERS Micro-SPE Clean-up Method for Pesticides in Feed Samples by GC-Orbitrap



PAL System μ -SPE Cartridges

Thank you very much for listening!

PAL SYSTEM
Ingenious sample handling

- Visit the new CTC μ -SPE cartridges at:
- <https://www.palsystem.com>
- <https://www.bgb-info.com/home.php?cat=38964>

PAL SYSTEM
Ingenious sample handling